

RABAC1 antibody

Catalog Number: 114491

Product name

RABAC1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RABAC1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in HEK-293 cells, human pancreas tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of RABAC1 recombinant protein (Accession Number: NM_006423). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114491(RABAC1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114491(RABAC1 antibody) at dilution of 1:50 (under 10x lens)