

RABGEF1 antibody



Catalog Number: 114498

Product name

RABGEF1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

RABGEF1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human brain tissue. Positive IF detected in PC-3 cells. Positive IHC detected in human endometrial cancer tissue. Observed molecular weight by Western-blot: 57kd




Preparation

This antibody was obtained by immunization of RABGEF1 recombinant protein (Accession Number: NM_001367735). Purification method: Antigen affinity purified.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114498(RABGEF1 antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:114498(RABGEF1 antibody) at dilution of 1:100 (under 10x lens)



RABGEF1 antibody



Catalog Number: 114498

Immunofluorescent analysis of PC-3 cells using
Catalog No:114498(RABGEF1 Antibody) at

dilution of 1:25 and Alexa Fluor 488-conjugated
AffiniPure Goat Anti-Rabbit IgG(H+L)