

# RABGEF1 antibody



Catalog Number: 114498

## Product name

RABGEF1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RABGEF1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human brain tissue. Positive IF detected in PC-3 cells. Positive IHC detected in human endometrial cancer tissue. Observed molecular weight by Western-blot: 57kd

## Validations



## Preparation

This antibody was obtained by immunization of RABGEF1 recombinant protein (Accession Number: NM\_001367735). Purification method: Antigen affinity purified.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114498(RABGEF1 antibody) at dilution of 1:300

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:114498(RABGEF1 antibody) at dilution of 1:100 (under 10x lens)



# RABGEF1 antibody



Catalog Number: 114498

---

Immunofluorescent analysis of PC-3 cells using  
Catalog No:114498(RABGEF1 Antibody) at

dilution of 1:25 and Alexa Fluor 488-conjugated  
AffiniPure Goat Anti-Rabbit IgG(H+L)