

RABGGTA antibody



Catalog Number: 114499

Product name

RABGGTA antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

RABGGTA Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells. Positive IHC detected in human lung tissue, human brain tissue. Observed molecular weight by Western-blot: 60 kDa




Preparation

This antibody was obtained by immunization of RABGGTA recombinant protein (Accession Number: NM_004581). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114499(RABGGTA antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

1:500


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human lung using Catalog No:114499(RABGGTA antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC

Immunohistochemical of paraffin-embedded human lung using Catalog No:114499(RABGGTA antibody) at dilution of 1:100 (under 40x lens)

Dilutions

Recommended Dilution:

WB: 1:500-1:5000