

RABL2A antibody



Catalog Number: 114500

Product name

RABL2A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RABL2A Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, human brain tissue. Positive IHC detected in human kidney. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of RABL2A recombinant protein (Accession Number: NM_001354423). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:114500(RABL2A antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114500(RABL2A antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114500(RABL2A antibody) at dilution of 1:50 (under 40x lens)