

RAC2 antibody

Catalog Number: 114506

Product name

RAC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAC2 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IP detected in HL-60 cells. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of RAC2 recombinant protein (Accession Number: NM_002872). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114506(RAC2 Antibody) at dilution of 1:600



IP Result of anti-RAC2 (IP:Catalog No:114506, 3ug; Detection:Catalog No:114506 1:300) with HL-60 cells lysate 4000ug.