

RACK1 antibody



Catalog Number: 114508

Product name

RACK1 antibody

Specificity

Human; other species not tested.

Antibody description

RACK1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human lung cancer tissue.

Preparation

This antibody was obtained by immunization of RACK1 recombinant protein (Accession Number: NM_006098). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:50-1:500

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114508(RACK1; GNB2L1 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:114508(RACK1; GNB2L1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)