

RAD23A antibody



Catalog Number: 114512

Product name

RAD23A antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

RAD23A Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 55kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of RAD23A recombinant protein (Accession Number: NM_005053). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114512(RAD23A antibody) at dilution of 1:500

Clonality