

# RAD51 antibody

Catalog Number: 114515

## Product name

RAD51 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RAD51 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, HeLa cells. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of RAD51 recombinant protein (Accession Number: BC001459). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114515(RAD51 Antibody) at dilution of 1:1000



IP Result of anti-RAD51 (IP:Catalog No:114515, 3ug; Detection:Catalog No:114515 1:1000) with mouse testis tissue lysate 4000ug.