

Rad51D antibody

Catalog Number: 114517

Product name

Rad51D antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Rad51D Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in COLO 320 cells, Jurkat cells. Observed molecular weight by Western-blot: 23-25kd

Preparation

This antibody was obtained by immunization of Rad51D recombinant protein (Accession Number: NM_002878). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:114517(Rad51D antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114517(Rad51D antibody) at dilution of 1:100 (under 10x lens)