

RAD9A antibody



Catalog Number: 114519

Product name

RAD9A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAD9A Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells, A549 cells, MCF7 cells. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of RAD9A recombinant protein (Accession Number: NM_004584). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114519(RAD9A antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:114519(RAD9A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)