

RAF1 antibody

Catalog Number: 114524

Product name

RAF1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAF1 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human prostate cancer tissue.

Preparation

This antibody was obtained by immunization of RAF1 recombinant protein (Accession Number: NM_002880). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114524(RAF1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114524(RAF1 Antibody) at dilution of 1:50 (under 40x lens)