

RALBP1 antibody

Catalog Number: 114531

Product name

RALBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RALBP1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HEK-293 cells, HepG2 cells, MCF7 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of RALBP1 recombinant protein (Accession Number: XM_047437281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114531(RALBP1 antibody) at dilution of 1:800



IP Result of anti-RALBP1 (IP:Catalog No:114531, 4ug; Detection:Catalog No:114531 1:800) with K-562 cells lysate 3600ug.