

# RALYL antibody



Catalog Number: 114532

## Product name

RALYL antibody

## Immunogen

[Human RALYL Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RALYL Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of RALYL recombinant protein (Accession Number: XM\_017013065). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114532(RALYL antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114532(RALYL Antibody) at dilution of 1:50 (under 40x lens)