

RAP2B antibody



Catalog Number: 114550

Product name

RAP2B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RAP2B Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 19-22 kDa

Preparation

This antibody was obtained by immunization of RAP2B recombinant protein (Accession Number: NM_002886). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114550(RAP2B antibody) at dilution of 1:500