

RAPGEF5 antibody



Catalog Number: 114551

Product name

RAPGEF5 antibody

Specificity

human, mouse; other species not tested.

Antibody description

RAPGEF5 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 68-75 kDa

Preparation

This antibody was obtained by immunization of RAPGEF5 recombinant protein (Accession Number: BC039203). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114551(RAPGEF5 antibody) at dilution of 1:500