

RCL1 antibody



Catalog Number: 114578

Product name

RCL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RCL1 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, A431 cells, HepG2 cells, MCF7 cells. Positive IF detected in A375 cells. Observed molecular weight by Western-blot: 41kd

Preparation

This antibody was obtained by immunization of RCL1 recombinant protein (Accession Number: NM_005772). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:114578(RCL1 antibody) at dilution of 1:500



Immunofluorescent analysis of A375 cells using Catalog No:114578(RCL1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)