

RCOR2 antibody



Catalog Number: 114580

Product name

RCOR2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RCOR2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells. Observed molecular weight by Western-blot: 58

Preparation

This antibody was obtained by immunization of RCOR2 recombinant protein (Accession Number: NM_173587). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114580(RCOR2 Antibody) at dilution of 1:1000