

# RDH11 antibody

Catalog Number: 114584

## Product name

RDH11 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RDH11 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of RDH11 recombinant protein (Accession Number: NM\_016026). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, IHC


## Dilutions

Recommended Dilution:


IHC: 1:50-1:500

IF: 1:10-1:100

## Validations



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:114584(RDH11 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunofluorescent analysis of PC-3 cells using Catalog No:114584(RDH11 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)