

# RBBP6 antibody

Catalog Number: 114592

## Product name

RBBP6 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RBBP6 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse placenta tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 200kd

## Preparation

This antibody was obtained by immunization of RBBP6 recombinant protein (Accession Number: NM\_032626). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114592(RBBP6 antibody) at dilution of 1:400



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114592(RBBP6 Antibody) at dilution of 1:500