

RBM15B antibody



Catalog Number: 114600

Product name

RBM15B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RBM15B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue. Observed molecular weight by Western-blot: 100-120 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of RBM15B recombinant protein (Accession Number: BC001367). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114600(RBM15B antibody) at dilution of 1:2000



IP Result of anti-RBM15B (IP:Catalog No:114600, 4ug; Detection:Catalog No:114600 1:500) with HeLa cells lysate 2500ug.