

RBM17 antibody



Catalog Number: 114603

Product name

RBM17 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RBM17 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells, MCF7 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of RBM17 recombinant protein (Accession Number: NM_032905). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114603(RBM17 antibody) at dilution of 1:800



Immunofluorescent analysis of HeLa cells, using RBM17 antibody Catalog No:114603 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).