

RBM28 antibody



Catalog Number: 114609

Product name

RBM28 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RBM28 Rabbit Polyclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 86 kDa

Preparation

This antibody was obtained by immunization of RBM28 recombinant protein (Accession Number: NM_018077). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114609(RBM28 antibody) at dilution of 1:1500



Immunofluorescent analysis of (10% Formaldehyde) fixed NIH/3T3 cells using Catalog No:114609(RBM28 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)