

# RBM4B antibody

Catalog Number: 114617

## Product name

RBM4B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RBM4B Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 36-40 kDa

## Preparation

This antibody was obtained by immunization of RBM4B recombinant protein (Accession Number: XM\_017018402). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114617(RBM4B antibody) at dilution of 1:500



IP Result of anti-RBM4B (IP:Catalog No:114617, 3ug; Detection:Catalog No:114617 1:300) with HEK-293 cells lysate 1000ug.