

RFC5 antibody



Catalog Number: 114626

Product name
RFC5 antibody

Specificity
Human,Mouse,Rat; other species not tested.

Antibody description
RFC5 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation
This antibody was obtained by immunization of RFC5 recombinant protein (Accession Number: NM_007370). Purification method: Antigen affinity purified.

Formulation
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage
Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type
Rabbit IgG

Applications
ELISA, WB

Dilutions
Recommended Dilution:
WB: 1:500-1:5000

Validations

-	+	-	+	Overexpression of RFC5-his-myc
+	+	-	-	Rabbit Anti human RFC5 polyclonal antibody
-	-	+	+	Mouse Anti his-tag monoclonal antibody

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114626(RFC5 Antibody) at dilution of 1:1000