

# RFTN2 antibody



Catalog Number: 114632

## Product name

RFTN2 antibody

IHC: 1:50-1:400


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

RFTN2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse brain tissue, mouse testis tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 56-60 kDa




## Preparation

This antibody was obtained by immunization of RFTN2 recombinant protein (Accession Number: NM\_144629). Purification method: Antigen affinity purified.

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114632(RFTN2 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:114632(RFTN2 Antibody) at dilution of 1:100 (under 10x lens).

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:114632(RFTN2 Antibody) at dilution of 1:100

# RFTN2 antibody



Catalog Number: 114632

---

(under 40x lens).