

RFX5 antibody

Catalog Number: 114635

Product name

RFX5 antibody

Specificity

Human; other species not tested.

Antibody description

RFX5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of RFX5 recombinant protein (Accession Number: XM_011509847). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114635(RFX5 antibody) at dilution of 1:400