

# RGS11 antibody



Catalog Number: 114641

## Product name

Polyclonal

RGS11 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Rat, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

RGS11 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, rat brain tissue. Observed molecular weight by Western-blot: 53 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of RGS11 recombinant protein (Accession Number: NM\_183337). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114641(RGS11 Antibody) at dilution of 1:1000