

RGS9 antibody



Catalog Number: 114650

Product name

RGS9 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RGS9 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 58kd, 70-75kd

Preparation

This antibody was obtained by immunization of RGS9 recombinant protein (Accession Number: BC022504). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114650(RGS9 antibody) at dilution of 1:500