

RHOC antibody



Catalog Number: 114660

Product name

RHOC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RHOC Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, C6 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of RHOC recombinant protein (Accession Number: NM_001042679). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114660(RHoc



antibody) at dilution of 1:500

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114660(RHoc antibody) at dilution of 1:500