

# RDH12 antibody



Catalog Number: 114668

## Product name

RDH12 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RDH12 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, COLO 320 cells. Observed molecular weight by Western-blot: 35kd

## Preparation

This antibody was obtained by immunization of RDH12 recombinant protein (Accession Number: NM\_152443). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114668(RDH12 antibody) at dilution of 1:200