

RDX antibody



Catalog Number: 114672

Product name

RDX antibody

Immunogen

[Human RDX Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RDX Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of RDX recombinant protein (Accession Number: XM_047427390). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114672(RDX Antibody) at dilution of 1:1000