

RECK antibody



Catalog Number: 114674

Product name

RECK antibody

Specificity

Human; other species not tested.

Antibody description

RECK Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, MCF7 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of RECK recombinant protein (Accession Number: NM_001316348). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114674(RECK antibody) at dilution of 1:600