

# REEP2 antibody



Catalog Number: 114680

## Product name

REEP2 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

human, mouse, rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

REEP2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 28kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of REEP2 recombinant protein (Accession Number: NM\_016606). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114680(REEP2 antibody) at dilution of 1:200

## Clonality