

RELB antibody

Catalog Number: 114684

Product name

RELB antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

IF: 1:10-1:100

Antibody description

RELB Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in Raji cells. Positive IF detected in HeLa cells. Positive IHC detected in human lymphoma tissue, human colon cancer tissue. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of RELB recombinant protein (Accession Number: NM_006509). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:114684(RELB Antibody) at dilution of 1:1000




IP Result of anti-RELB (IP:Catalog No:114684, 3ug; Detection:Catalog No:114684 1:300) with Raji cells lysate 2500ug.



RELB antibody

Catalog Number: 114684

Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:114684(RELB Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:114684(RELB Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG