

RELT antibody

Catalog Number: 114686

Product name

RELT antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RELT Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 46-50kd

Preparation

This antibody was obtained by immunization of RELT recombinant protein (Accession Number: NM_152222). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114686(RELT antibody) at dilution of 1:1000