

REM2 antibody



Catalog Number: 114687

Product name

REM2 antibody

Immunogen

[Human REM2 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

REM2 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of REM2 recombinant protein (Accession Number: BC035663). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114687(REM2 Antibody) at dilution of 1:1000