

RIBC2 antibody

Catalog Number: 114695

Product name

RIBC2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

human, mouse, rat; other species not tested.

Validations

Antibody description

RIBC2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human colon cancer. Positive IF detected in Hela. Observed molecular weight by Western-blot: 37kd




Preparation

This antibody was obtained by immunization of RIBC2 recombinant protein (Accession Number: XM_005261524). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114695(RIBC2 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114695(RIBC2 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded

RIBC2 antibody

Catalog Number: 114695

human colon cancer using Catalog No:114695(RIBC2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using RIBC2 antibody Catalog No:114695 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).