

# RIBC2 antibody



Catalog Number: 114695

## Product name

RIBC2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RIBC2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human colon cancer. Positive IF detected in Hela. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of RIBC2 recombinant protein (Accession Number: XM\_005261524). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114695(RIBC2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114695(RIBC2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded


# RIBC2 antibody



Catalog Number: 114695

---

human colon cancer using Catalog  
No:114695(RIBC2 antibody) at dilution of 1:50  
(under 40x lens)



Immunofluorescent analysis of Hela cells, using  
RIBC2 antibody Catalog No:114695 at 1:25  
dilution and Rhodamine-labeled goat anti-rabbit  
IgG (red).