

RIMKLB antibody



Catalog Number: 114705

Product name

RIMKLB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RIMKLB Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 43-48 kDa

Preparation

This antibody was obtained by immunization of RIMKLB recombinant protein (Accession Number: XM_017019695). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:114705(RIMKLB Antibody) at dilution of 1:300