

RIMKLB antibody



Catalog Number: 114705

Product name

Polyclonal

RIMKLB antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

RIMKLB Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 43-48 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of RIMKLB recombinant protein (Accession Number: XM_017019695). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:114705(RIMKLB Antibody) at dilution of 1:300