

RIPK2 antibody



Catalog Number: 114716

Product name

Polyclonal

RIPK2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

RIPK2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 46-55kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of RIPK2 recombinant protein (Accession Number: NM_003821). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114716(RIPK2 antibody) at dilution of 1:600