

# RKIP antibody



Catalog Number: 114718

## Product name

RKIP antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RKIP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse liver tissue, mouse testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of RKIP recombinant protein (Accession Number: NM\_002567). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:10-1:100

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114718(PEBP1 antibody) at dilution of 1:10



(under 10x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114718(PEBP1 Antibody) at dilution of 1:600