

RNASE11 antibody



Catalog Number: 114726

Product name

RNASE11 antibody

Specificity

Human; other species not tested.

Antibody description

RNASE11 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of RNASE11 recombinant protein (Accession Number: NM_001394191). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114726(RNASE11 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:114726(RNASE11 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:114726(RNASE11 antibody) at dilution of 1:50 (under 40x lens)