

RNASE13 antibody



Catalog Number: 114727

Product name

RNASE13 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNASE13 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse testis tissue. Observed molecular weight by Western-blot: 18-20 kDa

Preparation

This antibody was obtained by immunization of RNASE13 recombinant protein (Accession Number: NM_001012264). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114727(RNASE13 antibody) at dilution of 1:1500