

RNASEH2B antibody



Catalog Number: 114730

Product name

Polyclonal

RNASEH2B antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

RNASEH2B Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HeLa cells, MCF7 cells, mouse lung tissue. Observed molecular weight by Western-blot: 35kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of RNASEH2B recombinant protein (Accession Number: BC036744). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114730(RNASEH2B antibody) at dilution of 1:500