

# RNASEL antibody



Catalog Number: 114732

## Product name

RNASEL antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

Human, Rat, Mouse; other species not tested.

## Antibody description

RNASEL Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Observed molecular weight by Western-blot: 84 kDa

## Preparation

This antibody was obtained by immunization of RNASEL recombinant protein (Accession Number: NM\_021133). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IHC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114732(RNASEL Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:114732(RNASEL Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded


# RNASEL antibody



Catalog Number: 114732

---

human spleen tissue slide using Catalog No:114732(RNASEL Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:114732(RNASEL Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG