

RNF150 antibody



Catalog Number: 114742

Product name

RNF150 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RNF150 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human heart tissue, mouse brain tissue, mouse eye tissue, Y79 cells. Observed molecular weight by Western-blot: 43-44kd

Preparation

This antibody was obtained by immunization of RNF150 recombinant protein (Accession Number: BC101992). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114742(RNF150 antibody) at dilution of 1:500