

# RNF150 antibody



Catalog Number: 114742

## Product name

RNF150 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RNF150 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human heart tissue, mouse brain tissue, mouse eye tissue, Y79 cells. Observed molecular weight by Western-blot: 43-44kd

## Preparation

This antibody was obtained by immunization of RNF150 recombinant protein (Accession Number: BC101992). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114742(RNF150 antibody) at dilution of 1:500