

RNF167 antibody



Catalog Number: 114745

Product name

RNF167 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RNF167 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 46-50 kDa

Preparation

This antibody was obtained by immunization of RNF167 recombinant protein (Accession Number: NM_001370305). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:114745(RNF167 Antibody) at dilution of 1:300